



EXCEPTIONAL ENDOSCOPIC SOLUTIONS.

PRODUCT CATALOG

NEVER
COMPROMISE
ON QUALITY

DRAMATICALLY
IMPROVE YOUR
BOTTOM LINE

THROW AWAY
CONTRACTS
FOREVER

RELIABLE
SUPPLY
PARTNERSHIP

YOUR PARTNER
FOR GI ENDOSCOPY
DEVICES AND ACCESSORIES

A pair of blue nitrile gloves is shown against a light background. A white band is wrapped around the wrist of the glove in the foreground. The band has the text "MICRO-TECH" printed on it in a sans-serif font. The gloves have a textured surface with small ridges for grip. The background is a plain, light-colored surface.

MICRO-TECH

TABLE OF CONTENTS

BIOPSY	6
DILATION	10
POLYPECTOMY	15
HEMOSTASIS	20
GI ACCESSORIES	24
TERMS AND POLICIES	26
PRODUCT REFERENCE CHART	27
GENERAL TERMS AND CONDITIONS	30

MICRO-TECH ENDOSCOPY IS PROUD TO OFFER YOU A COMPLETELY UNIQUE, COMPLETELY EXCEPTIONAL WAY OF DOING BUSINESS.

Imagine working with an organization that puts all its resources into focusing on you. A company that gives you the freedom to choose on your own terms, and will never lock you into a contract.

We operate by a higher standard. We know that you're busy, value efficiency, and require perfection from the devices you use to keep your patients healthy.

At Micro-Tech, you have the freedom to purchase what you need, when you need it. The highest quality will never be compromised, and your satisfaction is always guaranteed.

With facilities in Asia, Europe, and America, we leverage our global reach to continually refine and provide the best products for endoscopic care.

You provide exceptional care. We provide exceptional support. Why settle for anything less?



Micro-Tech Endoscopy has teamed up with Henry Schein Medical in an exclusive, innovative partnership to deliver more value to our customers. You can now combine your Micro-Tech Endoscopy and Henry Schein orders for a true one-stop shopping experience.

CURRENTLY SERVING OVER
3,500 ASCs

NEXT DAY DELIVERY*
LATE DAY ORDERING CUT-OFF TIMES
NOMINAL SHIPPING AND HANDLING CHARGES**
ORDER ITEMS IN LOW UNITS OF MEASURE†

HOW TO ORDER



CALL

877.552.4027
Monday – Friday
8:30 AM – 5:30 PM EST



EMAIL

customer care@mtendoscopy.com



FAX

734.272.4160



MAIL

Micro-Tech Endoscopy USA
2855 Boardwalk Drive
Ann Arbor, MI 48104

*For majority of customers – some restrictions apply

**Minimum purchase of \$200.00

†Excludes bite blocks

BIOPSY

BIOPSY FORCEPS

With our wide selection of cost-effective disposable biopsy forceps, you are perfectly equipped for every procedure.

- Available jaw diameters of 1.8mm, 2.3mm and 3.0mm, and are offered in Oval, Serrated, Alligator jaws, with or without a spike
- Available lengths of 160cm, 180cm and 230mm.
- Standard Oval and Serrated offer larger capacity
- All forceps feature entry and exit indicator markings.

KEY BENEFITS

ERGONOMIC HANDLE

Three-ring comfort with cushioned lining and short throw provides improved grip and accommodates smaller hands.

DIRECT DRIVE SYSTEM

The single drive wire and hinge linkage delivers smooth activation around curves and tortuous anatomy. The system delivers more closure force to the jaws which provides for improved cutting ability and assures complete jaw closure.

SCOPE FRIENDLY DESIGN

The design of the Micro-Tech forceps family is based upon the robust and proven design of original reusable biopsy instruments. These designs provide optimal bite, while minimizing the risk of scope damage.

BIOPSY FORCEP FOR 2.0MM SCOPE CHANNELS



Order Number	Henry Schein Item Number	Description	Jaw Style	Sheath Color	Working Length (cm)	Minimum Channel Size (mm)	Package Units
BF40021	129-0175	Mini Forceps	Oval	Yellow	160	2	10/Box
BF40051	129-0176	Mini Forceps with Spike	Oval	Yellow	160	2	10/Box
BF40141	129-0181	Mini Forceps	Alligator	Yellow	160	2	10/Box

BIOPSY FORCEP FOR 2.8MM SCOPE CHANNELS

STERILE

②

Order Number	Henry Schein Item Number	Description	Jaw Style	Sheath Color	Working Length (cm)	Minimum Channel Size (mm)	Package Units
BF40081	129-0179	Standard Forceps	Oval	Orange	180	2.8	10/Box
BF40091	129-0677	Standard Forceps	Oval	Blue	230	2.8	10/Box
BF40096	128-5668	Standard Forceps	Oval	Blue	230	2.8	50/Box
BF40111	129-0180	Standard Forceps with Spike	Oval	Orange	180	2.8	10/Box
BF40121	129-0678	Standard Forceps with Spike	Oval	Blue	230	2.8	10/Box
BF40126	128-5669	Standard Forceps with Spike	Oval	Blue	230	2.8	50/Box
BF40191	129-0184	Standard Forceps	Alligator	Orange	180	2.8	10/Box
BF40201	129-0679	Standard Forceps	Alligator	Blue	230	2.8	10/Box
BF40206	128-5670	Standard Forceps	Alligator	Blue	230	2.8	50/Box
BF40221	129-0185	Standard Forceps with Spike	Alligator	Orange	180	2.8	10/Box
BF40231	129-0680	Standard Forceps with Spike	Alligator	Blue	230	2.8	10/Box
BF40236	128-5671	Standard Forceps with Spike	Alligator	Blue	230	2.8	50/Box
BF40431	132-9020	Standard Forceps	Serrated	Orange	180	2.8	10/Box
BF40441	132-9017	Standard Forceps	Serrated	Blue	230	2.8	10/Box
BF40446	131-2212	Standard Forceps	Serrated	Blue	230	2.8	50/Box
BF40411	132-9018	Standard Forceps with Spike	Serrated	Orange	180	2.8	10/Box
BF40421	132-9019	Standard Forceps with Spike	Serrated	Blue	230	2.8	10/Box
BF40426	131-2211	Standard Forceps with Spike	Serrated	Blue	230	2.8	50/Box

BIOPSY FORCEP FOR 3.2MM SCOPE CHANNELS

STERILE

②

Order Number	Henry Schein Item Number	Description	Jaw Style	Sheath Color	Working Length (cm)	Minimum Channel Size (mm)	Package Units
BF40061	129-0177	Jumbo Biopsy	Oval	Blue	230	3.2	10/Box
BF40071	129-0178	Jumbo Biopsy with Spike	Oval	Blue	230	3.2	10/Box
BF40171	129-0182	Jumbo Biopsy	Alligator	Blue	230	3.2	10/Box
BF40181	129-0183	Jumbo Biopsy with Spike	Alligator	Blue	230	3.2	10/Box

BIOPSY FORCEP HOT

STERILE

②

Order Number	Henry Schein Item Number	Description	Jaw Style	Sheath Color	Working Length (cm)	Minimum Channel Size (mm)	Package Units
HF48521	129-0681	Biopsy Forceps Hot Lower	Oval	Blue	230	2.8	10/Box
HF48541	131-2213	Biopsy Forceps Hot Lower	Serrated	Blue	230	2.8	10/Box



LARGE CAPACITY BIOPSY FORCEPS

Twin-wire mechanism and precision stamping technology for assured tissue sampling

- Pivoting jaws to aid in sampling where space may be limited
- Elongated Jaw for greater tissue retention
- With needle option for multiples bites in a single pass
- Entry - Exit indicator markings on sheath
- Flexible and pushable smooth catheter

BIOPSY FORCEP FOR 2.8MM SCOPE CHANNELS

STERILE 2

Order Number	Henry Schein Item Number	Description	Jaw Style	Sheath Color	Working Length (cm)	Minimum Channel Size (mm)	Package Units
BF40286	133-4667	Large Capacity Forceps	Alligator	Blue	230	2.8	50/Box
BF40306	133-4666	Large Capacity Forceps with Spike	Alligator	Blue	230	2.8	50/Box



GI CYTOLOGY BRUSHES

Stiff brushes for reliable tissue sampling

CYTOLOGY BRUSH

STERILE 2

Order Number	Henry Schein Item Number	Description	Brush O.D. (mm)	Catheter O.D. (mm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
CY49021	133-5541	Cytology Brush	2.0	1.8	180	2.0	20/Box
CY49051	133-5542	Cytology Brush	3.0	2.3	240	2.8	20/Box

PRECISION & CONTROL

RETRIEVAL FORCEPS

Single-use grasping forceps for secure removal of tissue or foreign bodies.

- 3 Jaw types
- Ergonomic 3-ring handle
- Entry - Exit markers
- Sterile - Single Use Convenience

RETRIEVAL FORCEPS

STERILE 2

Order Number	Henry Schein Item Number	Description	Jaw Opening (mm)	Color	Working Length (mm)	Minimum Channel Size (mm)	Package Units
GF61021	131-2207	Rat Tooth	6.3	Blue	230	2.8	10/Box
GF61041	131-2208	Alligator Tooth	8.1	Blue	230	2.8	10/Box
GF61061	131-2209	Mixed Tooth	8.1	Blue	230	2.8	10/Box

DILATION

STRENGTH AND CONTROL MULTI-STAGE DILATION BALLOONS

Choose strength and control when dilating strictures with a dilation balloon.

- Multi-Stage balloons feature a premium clear Pebax® balloon with rounded shoulders for balloon endoscopy
- Fixed wire esophageal balloons offer fast, reliable performance
- Pre-loaded .035 wire-guided balloons assist placement in tortuous or altered anatomy

KEY BENEFITS

CHOICE OF MATERIAL

The high quality Pebax material produces a balloon with the features and performance your providers demand. Accurate three-size radial expansion, rounded shoulders, clear visibility and strength with control.

VISIBILITY UNDER FLUOROSCOPY

All wire-guided balloons feature two radiopaque markers at the balloon shoulder, allowing for accurate placement across the stricture. The dedicated inflation port permits the introduction of contrast to further enhance visibility if desired.

TAPERED TIP

The rounded tip aids in the smooth passage of the balloon across the stricture. The wire-guided balloon is provided with a pre-loaded .035 guidewire that has a soft, flexible tip to assist in placement in tortuous or altered anatomy.

WIRE-GUIDED ESOPHAGEAL, PYLORIC, COLONIC MULTI-STAGE BALLOONS

STERILE

2

Order Number	Henry Schein Item Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
MB23031	125-2734	Reliant™ Multi-Stage Dilatation Balloon	6-7-8	5.5	230	2.8	1/Box
MB23071	125-2736	Reliant™ Multi-Stage Dilatation Balloon	8-9-10	5.5	230	2.8	1/Box
MB23111	125-2738	Reliant™ Multi-Stage Dilatation Balloon	10-11-12	5.5	230	2.8	1/Box
MB23151	125-2740	Reliant™ Multi-Stage Dilatation Balloon	12-13.5-15	5.5	230	2.8	1/Box
MB23191	125-2742	Reliant™ Multi-Stage Dilatation Balloon	15-16.5-18	5.5	230	2.8	1/Box
MB23221	125-2744	Reliant™ Multi-Stage Dilatation Balloon	18-19-20	5.5	230	2.8	1/Box

FIXED WIRE ESOPHAGEAL MULTI-STAGE BALLOONS

STERILE

2

Order Number	Henry Schein Item Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24015	129-0163	Rapide™ Esophageal Dilatation Balloon	6-7-8	8	180	2.8	2/Box
FW24025	129-0164	Rapide™ Esophageal Dilatation Balloon	8-9-10	8	180	2.8	2/Box
FW24035	129-0165	Rapide™ Esophageal Dilatation Balloon	10-11-12	8	180	2.8	2/Box
FW24045	129-0166	Rapide™ Esophageal Dilatation Balloon	12-13.5-15	8	180	2.8	2/Box
FW24055	129-0167	Rapide™ Esophageal Dilatation Balloon	15-16.5-18	8	180	2.8	2/Box
FW24065	129-0168	Rapide™ Esophageal Dilatation Balloon	18-19-20	8	180	2.8	2/Box

INFLATION DEVICE

STERILE

2

Order Number	Henry Schein Item Number	Description	Package Units
DB25011	131-4493	Micro-Tech Balloon Inflation Device 60 ml with lock	1/Box
QL6015	125-2750	Atrion 60cc, Single Use	1/Box



ABC™ 3-STEP INFLATION BALLOONS

- Offers the radial strength to accurately and consistently dilate strictures
- Easily inflated in three steps, eliminating the need for multiple balloons for the same procedure
- An economical solution vs. typical radial expansion balloons
- Available as a flexible fixed wire esophageal 8cm balloon on a 6Fr catheter, featuring rapid inflation/deflation performance, or as a wire-guided esophageal, pyloric, colonic 5.5cm balloon, on a dual channel 7Fr catheter with a .035" flexible tip guide-wire

KEY BENEFITS

CHOICE OF MATERIAL

The durable nylon material delivers strength with accurate sizing. The non-compliant material provides the radial force required for consistent dilation. Inflated in 3 steps, the ABC balloon provides an affordable alternative to radial expanding balloons.

VISIBILITY UNDER FLUOROSCOPY

All wire guided balloons feature two radiopaque markers at the balloon shoulder, allowing for accurate placement across the stricture. The dedicated inflation port permits the introduction of contrast to further enhance visibility if desired.

TAPERED TIP

The rounded tip aids in the smooth passage of the balloon across the stricture. The wire guided balloon is provided with a pre-loaded .035 guide-wire that has a soft, flexible tip to assist in placement in tortuous or altered anatomy.

WIRE-GUIDED ESOPHAGEAL, PYLORIC, COLONIC BALLOONS

STERILE

②

Order Number	Henry Schein® Item Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
BD20091	125-2689	TTS - Wire Guided Dilation Balloon	6	5.5	240	2.8	1/Box
BD20181	125-2691	TTS - Wire Guided Dilation Balloon	8	5.5	240	2.8	1/Box
BD20271	125-2693	TTS - Wire Guided Dilation Balloon	10	5.5	240	2.8	1/Box
BD20361	125-2696	TTS - Wire Guided Dilation Balloon	12	5.5	240	2.8	1/Box
BD20451	125-2698	TTS - Wire Guided Dilation Balloon	14	5.5	240	2.8	1/Box
BD20631	125-2700	TTS - Wire Guided Dilation Balloon	16	5.5	240	2.8	1/Box
BD20721	125-2702	TTS - Wire Guided Dilation Balloon	18	5.5	240	2.8	1/Box

FIXED WIRE ESOPHAGEAL 3-STEP INFLATION BALLOONS

STERILE

②

Order Number	Henry Schein Item Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24075	128-5660	ABC™ Esophageal Dilation Balloon	6	8	180	2.8	2/Box
FW24085	128-5661	ABC™ Esophageal Dilation Balloon	8	8	180	2.8	2/Box
FW24095	128-5662	ABC™ Esophageal Dilation Balloon	10	8	180	2.8	2/Box
FW24105	128-5663	ABC™ Esophageal Dilation Balloon	12	8	180	2.8	2/Box
FW24125	128-5664	ABC™ Esophageal Dilation Balloon	15	8	180	2.8	2/Box
FW24145	128-5665	ABC™ Esophageal Dilation Balloon	18	8	180	2.8	2/Box
FW24155	128-5666	ABC™ Esophageal Dilation Balloon	20	8	180	2.8	2/Box

INFLATION DEVICE

STERILE

②

Order Number	Henry Schein Item Number	Description	Package Units
DB25011	131-4493	Micro-Tech Balloon Inflation Device 60 ml with lock	1/Box
QL6015	125-2750	Atrion 60cc, Single Use	1/Box

SIMPLE SOLUTION BALLOON INFLATION DEVICE

The Micro-Tech Inflation Device was designed to accommodate smaller hands, featuring a lever-lock.



- Lever-lock controlled with the holding hand
- Plunger moved to prime, inflate, or create a vacuum
- Luminous dial includes atmosphere on the outer scale for clear visualization
- Priced competitively
- 60mL single-use syringe includes a fill mark for Micro-Tech Balloons, in combination delivers significant value



LOCKED



UNLOCKED

INFLATION DEVICE

STERILE 2


Order Number	Henry Schein Item Number	Description	Package Units
DB25011	131-4493	Micro-Tech Balloon Inflation Device 60 ml with lock	1/Box
QL6015	125-2750	Atrion 60cc, Single Use	1/Box

POLYPECTOMY

PRECISE CUT HOT SNARES

Our line of hot and cold snares features effective cutting wires and user-friendly handles at cost-effective prices.

- With the innovative power connector, the snare can be operated and rotated with ease.
- Our hot snares are available in a compliant braided wire or a firmer wire for effective tissue acquisition and resection



KEY BENEFITS

SHORTER THROW HANDLE

The handle throw is adjusted according to the snare size. Small and Medium snares have a shorter throw than Large or Extra Large snares. A shorter throw improves comfort for users with smaller hands.

BRAIDED WIRE

The braided wire combines the tissue acquisition features of a firm snare with the flexibility and smooth deployment of a soft snare. Suitable for daily use in facilities looking to minimize inventory.

INNOVATIVE POWER CONNECTION

The handle features a dedicated power connection, separate from the finger ring slider. This feature permits unrestricted movement and rotation of the snare handle and may reduce procedural delays.

HOT SNARES BRAIDED



Order Number	Henry Schein Item Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
PS51011	129-0186	Mini Braided Oval	10	7	230	2.8	10/Box
PS51021	129-0673	Small Braided Oval	15	7	230	2.8	10/Box
PS51031	129-0674	Medium Braided Oval	20	7	230	2.8	10/Box
PS51041	129-0675	Large Braided Oval	24	7	230	2.8	10/Box
PS51051	129-0187	Extra Large Braided Oval	30	7	230	2.8	10/Box

HOT SNARES FIRM



Order Number	Henry Schein Item Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
PS51071	129-0188	Mini Oval Firm	10	7	230	2.8	10/Box
PS51081	129-0189	Small Oval Firm	15	7	230	2.8	10/Box
PS51091	129-0190	Medium Oval Firm	20	7	230	2.8	10/Box
PS51101	129-0191	Large Oval Firm	24	7	230	2.8	10/Box
PS51111	129-0192	Extra Large Oval Firm	30	7	230	2.8	10/Box
PS51121	129-0193	Jumbo Oval Firm	36	7	230	2.8	10/Box

HOT SNARES HEX FIRM



Order Number	Henry Schein Item Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
PS52031	129-0194	Small Hex Firm	15	7	230	2.8	10/Box
PS52041	129-0195	Large Hex Firm	25	7	230	2.8	10/Box

EFFICIENT RESECTION

COLD SNARES

- Designed specifically to resect diminutive polyps
- Available in two sizes – 10mm and 15mm
- Thin cutting wire provides for a guillotine-type cut
- Loop shape assists in capturing smaller polyps
- Stiffer catheter aids in insertion and provides support during resection

KEY BENEFITS

DEDICATED COLD SNARE

Specifically designed for resecting diminutive polyps using a cold snare technique, the oval or diamond shape loop and flexible wire, capture smaller polyps with ease.

THIN WIRE

The thin 0.23 mm cutting wire provides a guillotine cut that reduces the force needed to complete the resection and may limit tearing or snagging.

COMFORTABLE ERGONOMIC HANDLES

Accommodating various hand sizes, the ergonomic three-ring handle features a shortened throw and comfortable grip. The stiff catheter provides support for a smooth effortless cut.

COLD SNARES




Order Number	Henry Schein Item Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
CS50011	129-0671	Diamond Thin Wire	10	7	230	2.8	10/Box
CS50021	129-0672	Diamond Thin Wire	15	7	230	2.8	10/Box
CS50031	129-0160	Oval Thin Wire	10	7	230	2.8	10/Box
CS50041	129-0161	Oval Thin Wire	15	7	230	2.8	10/Box

TRUE BLUE BLUE BEACON™ G.I. MARKER



- 0.4% Indigo Carmine solution
- Useful in GI procedures that require a contrast stain to temporarily mark or highlight tissue
- Can be diluted and introduced to the irrigation, sprayed or injected around the area of interest
- Available in a convenient luer lock 5ml syringe. The 0.4% Indigo Carmine solution is supplied sterile, 5 syringes per box.



BLUE BEACON™ G.I. MARKER

STERILE



Order Number	Henry Schein Item Number	Description	Package Units
IC62021	128-5667	Blue Beacon GI Marker 5ml Syringe Indigo Carmine 0.4% Sterile	5/Box

This solution is NOT FOR INTRAVENOUS OR INTRAMUSCULAR INJECTION.
It is a single use device indicated for temporary marking of the GI tract.

SPRAY CATHETER

STERILE



Order Number	Henry Schein Item Number	Description	Working Length (cm)	Minimum Channel Size (mm)	Package Units
SC62241	133-5544	Spray Catheter	220	2.2	20/Box
SC62301	133-5545	Spray Catheter	250	2.8	20/Box

ENTRAP™ POLYP RETRIEVAL NET

The Entrap retrieval system delivers a quick and efficient retrieval solution for polyps and tissue for histopathologic review. Utilizing advanced materials and production techniques, the Entrap facilitates improved consistency and control of the net deployment, position, and closure.



KEY BENEFITS

REFINED NET

Improved visibility through the net material for identification of the retrieval target. Durable enough for repeated deployment, but soft to gently capture and preserve the specimen.

ADVANCED PROCESSES

Cinching, tearing, and binding of the net material is reduced by the minimal net selvage and bonded net attachment. This increases efficiency and minimizes the time required to reliably capture specimens.

SMOOTH AND STRONG

A stiffer catheter material improves device insertion into the endoscope and positioning of the net. The additional columnar strength allows for gentle movement of a specimen for optimal capture.

ENTRAP RETRIEVAL SYSTEM



Order Number	Henry Schein Item Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
NT54511	133-5543	Entrap Polyp Retrieval Net	30	7	230	2.8	10/Box


HEMOSTASIS

CLIP WITH CONFIDENCE **SURECLIP®**

Clips need to be reliable. They need to be accurate. They need to allow you the flexibility to reposition or rotate as much as is required to deliver better outcomes.


Accurate positioning prior to deployment can reduce both procedure time and the number of clips needed to achieve hemostasis. SureClip achieves this by design delivering outstanding repositionability and reliable rotation prior to deployment, in various scope positions. SureClip's short stem aids placement in narrow lumen and improves visibility.

The proprietary clip design provides reliable deployment, may improve retention, and offers a choice of jaw sizes.



SURECLIP

Shortest stem, approximately 5mm,
for use in narrow lumen



SURECLIP^{PLUS}

Widest jaw opening for confident defect closure



NEW PRODUCT!

SURECLIP^{MINI}

Lowest jaw profile for unobstructed view during placement



KEY BENEFITS

REPOSITIONING

SureClip's unique design permits opening and closing the jaw prior to deployment. Being able to reposition a clip may help improve placement accuracy. Fenestrations on the SureClip^{PLUS}, accommodate tissue and may enhance retention.

RELIABLE FULL ROTATION

SureClip can be rotated, helping to provide the correct orientation for tissue approximation or vessel closure. The rotation handle on the SureClip^{PLUS} improves performance and enhances the user experience.

SHORT STEM

A shorter stem makes the clip less obtrusive, improving visualization of the target area, particularly when multiple clips are placed in close proximity.



SURECLIP^{MINI} LOWEST PROFILE

STERILE

2

Order Number	Henry Schein Item Number	Opening Width (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size(mm)	Package Units
RC30415	132-5180	8	Max 2.6	235	2.8	2/Box
RC30411	132-5187	8	Max 2.6	235	2.8	10/Box



SURECLIP SHORTEST STEM

STERILE

2

RC30445	132-5724	11	Max 2.6	235	2.8	2/Box
RC30441	132-5723	11	Max 2.6	235	2.8	10/Box



SURECLIP^{PLUS} WIDEST JAW OPENING

STERILE

2

RC30385	128-5657	16	Max 2.6	235	2.8	2/Box
RC30381	128-5655	16	Max 2.6	235	2.8	10/Box

*Jaws may not open to maximum reach until placed in contact with tissue

SHARP INNOVATION INJECTION NEEDLES



KEY BENEFITS

ONE-HANDED OPERATION

Ergonomically designed one-handed operation enables confident needle deployment that enables rapid needle retraction.

PATENTED RATCHET DESIGN

Provides audible confirmation and secure locking of the deployed needle during use. As the needle is advanced, the ratchet gives tactile feedback and helps to reduce needle bounce.

PROTECTIVE METAL HUB

The rounded metal hub protects and guides the needle during insertion and deployment. The hub and firm catheter may support use for tamponade.

INJECTION NEEDLES



Order Number	Henry Schein Item Number	Description	Needles Gauge/Length (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
IN33131	129-0169	Injection Needle Upper	22/5	2.3	180	2.8	10/Box
IN33221	129-0171	Injection Needle Upper	25/5	2.3	180	2.8	10/Box
IN33151	129-0170	Injection Needle Lower	22/5	2.3	230	2.8	10/Box
IN33241	129-0172	Injection Needle Lower	25/5	2.3	230	2.8	10/Box

EFFECTIVE LIGATION HEMORRHOIDS MULTIPLE BAND LIGATOR SET

Endoscopic banding has been demonstrated as an effective and simple treatment for internal hemorrhoids. Banding is readily performed in an ambulatory setting and can result in less pain and have a shorter recovery period than surgical hemorrhoidectomy.

The Micro-Tech Hemorrhoids Multiple Band Ligator Set is a cost-effective solution for the outpatient treatment of hemorrhoids during a routine colonoscopy.



KEY BENEFITS

EXCELLENT VISIBILITY

A clear ligating barrel with rear mounted bands optimizes the endoscopic view while maintaining suction.

BAND DESIGN

The square band design provides consistent shape and reliable retention. The set comes with four pre-loaded latex-free bands.

USER-FRIENDLY SYSTEM

The handle unit tie hole enables quick and easy loading of the pull wire. The alert band notifies the operator when only one band remains.

HEMORRHOIDS MULTIPLE BAND LIGATOR SET

STERILE

2

Order Number	Henry Schein Item Number	Description	Minimum Working Channel (mm)	Working Length (cm)	Package Units
HL32161	131-9042	Hemorrhoids Multiple Band Ligator Set, 4 Latex Free Bands, compatible with 9.4 - 13 mm scope O.D.	2.8	220	1/Box

GI ACCESSORIES

PROTECT, COLLECT, CLEAN.

PROTECT



BITE BLOCKS

Micro-Tech Bite Blocks offer comfort and protection during upper endoscopic procedures. Available in 38Fr, 60Fr and 60Plus sizes, with or without dental rim.

Order Number	Henry Schein Item Number	Description	Package Units
BB60011	125-2687	Bite Block, Adult, With Strap	100/Box
BB60021	125-2688	Bite Block, Pediatric, With Strap	100/Box
BB60041	131-2210	Bite Block, Adult w/o dental rim. 60+Fr with strap	100/Box
BB60051	133-5546	Bite Block, Adult, 60+ Fr with strap	100/Box



CATCHALL SCOPE TRANSPORT PAD

The CatchAll Transport Pad offers compliance and protection while transporting contaminated scopes.

Order Number	Henry Schein Item Number	Description	Package Units
CP62121	131-2218	CatchAll Scope Transport Pad	25/Box

GI ACCESSORIES

COLLECT



POLYP TRAPS ²

Polyp Retrieval Traps offer improved efficiency, reduced handling and cost savings.

The unique Quad trap features removable collection trays and a secondary filter for peace of mind.

Order Number	Henry Schein Item Number	Description	Package Units
PC54011	129-0676	Polyp Trap with 4 Chambers	25/Box
PC54021	131-2219	Polyp Trap single chamber	25/Box
PC54031	131-2214	Quad Polyp Trap with removable chambers	25/Box

CLEAN



CLEANING BRUSHES ²

Micro-Tech's Cleaning Brushes combine performance and compliance. Available in combination kits, the brush range features 6/10mm nylon brush heads designed for effective gentle cleaning.

Order Number	Henry Schein Item Number	Description	Package Units
QL22011	131-2220	Single ended 6mm channel brush	50/Box
QL22021	131-2222	Double ended 6mm channel brush	50/Box
QL22031	131-2221	Single ended 6mm channel and attached 10mm valve brush	50/Box
QL22041	131-2215	Double ended 5mm-10mm valve brush	50/Box
QL22051	131-2216	Single ended 6mm channel brush + 5mm-10mm valve brush	50/Box
QL22061	131-2217	Double ended 6mm channel brush + 5mm-10mm valve brush	50/Box

TERMS AND POLICIES

SHIPMENT AND PAYMENT TERMS

Items are normally shipped Federal Express Ground or UPS Ground, FOB shipping location. Shipping charges will be added to the invoice. Payment Terms are net 30 days. All prices are subject to change without notice. Express deliveries of stock items can be arranged and will be subject to additional freight charges.

SPECIAL ORDER ITEMS

In addition to stock items listed within this catalog, Micro-Tech Endoscopy manufactures non-stock items. Non-stock items may be subject to longer lead times. Contact Customer Support prior to ordering for more information.

STANDING ORDER SHIPMENTS

A standing order shipment simplifies ordering and shipments of regularly purchased products. Products ordered routinely on a regular basis can be scheduled to ship at an agreed upon, regular interval. Standing order shipments offer convenience and can be canceled at any time.

RETURNED GOODS POLICY

Micro-Tech Endoscopy offers a 100% refund on all unused and unopened products, in the original product box, that are returned within 90 days from the date of purchase. Written authorization must be obtained from Micro-Tech Endoscopy prior to return.

PRODUCT STORAGE

Micro-Tech Endoscopy products should be stored in a cool, dry, clean, well ventilated place. Do not expose the packaging to organic solvents, ionizing radiation or ultraviolet radiation.

DESIGN CHANGES

Micro-Tech Endoscopy reserves the right to discontinue or change products without notice. If any item has been discontinued, alternatives will be offered.

UPN DESCRIPTION

Now products may be ordered using the UPN (Unique Part Number) as well as the traditional Ref Number.

QUALITY POLICY

Micro-Tech Endoscopy is committed to supplying a comprehensive portfolio of quality, cost-effective gastrointestinal devices. We are committed to maintaining a quality system that provides for safe and effective products and services that meet the needs of our patients, customers and stakeholders.

PRODUCT REFERENCE CHART

Order Number	Ref #	UPN Pouch	Henry Schein Item Number	Product Family
BB60011	AC01-103.A	BB60010	125-2687	GI Accessories
BB60021	AC01-103.P	BB60020	125-2688	GI Accessories
BB60041	AC01-108.M	BB60040	131-2210	GI Accessories
BB60051	AC01-109.M	BB60050	133-5546	GI Accessories
BD20091	BDC-6/55-7/24-A	BD20090	125-2689	Dilation
BD20181	BDC-8/55-7/24-A	BD20180	125-2691	Dilation
BD20271	BDC-10/55-7/24-A	BD20270	125-2693	Dilation
BD20361	BDC-12/55-7/24-A	BD20360	125-2696	Dilation
BD20451	BDC-14/55-7/24-A	BD20450	125-2698	Dilation
BD20631	BDC-16/55-7/24-A	BD20630	125-2700	Dilation
BD20721	BDC-18/55-7/24-A	BD20720	125-2702	Dilation
BF40021	NBF02-11018160	BF40020	129-0175	Biopsy
BF40051	NBF02-11118160	BF40050	129-0176	Biopsy
BF40061	NBF03-01030230	BF40060	129-0177	Biopsy
BF40071	NBF03-01130230	BF40070	129-0178	Biopsy
BF40081	NBF03-11023180	BF40080	129-0179	Biopsy
BF40091	NBF03-11023230	BF40090	129-0677	Biopsy
BF40096	NBF03-11023230	BF40094	128-5668	Biopsy
BF40111	NBF03-11123180	BF40110	129-0180	Biopsy
BF40121	NBF03-11123230	BF40120	129-0678	Biopsy
BF40126	NBF03-11123230	BF40124	128-5669	Biopsy
BF40141	NBF12-11018160	BF40140	129-0181	Biopsy
BF40171	NBF13-01030230	BF40170	129-0182	Biopsy
BF40181	NBF13-01130230	BF40180	129-0183	Biopsy
BF40191	NBF13-11023180	BF40190	129-0184	Biopsy
BF40201	NBF13-11023230	BF40200	129-0679	Biopsy
BF40206	NBF13-11023230	BF40204	128-5670	Biopsy
BF40221	NBF13-11123180	BF40220	129-0185	Biopsy
BF40231	NBF13-11123230	BF40230	129-0680	Biopsy
BF40236	NBF13-11123230	BF40234	128-5671	Biopsy
BF40286	LCNBF13-10024230	BF40284	133-4667	Biopsy
BF40306	LCNBF13-10124230	BF40304	133-4666	Biopsy
BF40411	NBF63-11123180	BF40410	132-9018	Biopsy
BF40421	NBF63-11123230	BF40420	132-9019	Biopsy
BF40426	NBF63-11123230	BF40424	131-2211	Biopsy
BF40431	NBF63-11023180	BF40430	132-9020	Biopsy
BF40441	NBF63-11023230	BF40440	132-9017	Biopsy
BF40446	NBF63-11023230	BF40444	131-2212	Biopsy
CP62121	CTP-01025	CP62120	131-2218	GI Accessories
CS50011	CS3-11023230	CS50010	129-0671	Polypectomy
CS50021	CS3-11523230	CS50020	129-0672	Polypectomy
CS50031	CS3-21023230	CS50030	129-0160	Polypectomy
CS50041	CS3-21523230	CS50040	129-0161	Polypectomy
CY49021	CYB-20S-18C	CY49020	133-5541	Biopsy
CY49051	CYB-30S-24C	CY49050	133-5542	Biopsy
DB25011	BID-15-60-01	DB25010	131-4493	Dilation

Order Number	Ref #	UPN Pouch	Henry Schein Item Number	Product Family
FW24015	FWB-06/80-6/18-M	FW24014	129-0163	Dilation
FW24025	FWB-08/80-6/18-M	FW24024	129-0164	Dilation
FW24035	FWB-10/80-6/18-M	FW24034	129-0165	Dilation
FW24045	FWB-12/80-6/18-M	FW24044	129-0166	Dilation
FW24055	FWB-15/80-6/18-M	FW24054	129-0167	Dilation
FW24065	FWB-18/80-6/18-M	FW24064	129-0168	Dilation
FW24075	FWB-06/80-6/18-S	FW24074	128-5660	Dilation
FW24085	FWB-08/80-6/18-S	FW24084	128-5661	Dilation
FW24095	FWB-10/80-6/18-S	FW24094	128-5662	Dilation
FW24105	FWB-12/80-6/18-S	FW24104	128-5663	Dilation
FW24125	FWB-15/80-6/18-S	FW24124	128-5664	Dilation
FW24145	FWB-18/80-6/18-S	FW24144	128-5665	Dilation
FW24155	FWB-20/80-6/18-S	FW24154	128-5666	Dilation
GF61021	NGF03-11023230	GF61020	131-2207	Biopsy
GF61041	NGF13-11023230	GF61040	131-2208	Biopsy
GF61061	NGF33-11023230	GF61060	131-2209	Biopsy
HF48521	HBf55-11023230	HF48520	129-0681	Biopsy
HF48541	HBf65-11023230	HF48540	131-2213	Biopsy
HL32161	HMBLS-4F-NL	HL32160	131-9042	Hemostasis
IC62021	P-DYZ-04-05	IC62020	128-5667	Polypectomy
IN33131	IN12-225231802	IN33130	129-0169	Hemostasis
IN33151	IN12-225232302	IN33150	129-0170	Hemostasis
IN33221	IN12-255231802	IN33220	129-0171	Hemostasis
IN33241	IN12-255232302	IN33240	129-0172	Hemostasis
MB23031	MBD-AM-BB-2	MB23030	125-2734	Dilation
MB23071	MBD-BM-BB-2	MB23070	125-2736	Dilation
MB23111	MBD-CM-BB-2	MB23110	125-2738	Dilation
MB23151	MBD-DM-BB-2	MB23150	125-2740	Dilation
MB23191	MBD-EM-BB-2	MB23190	125-2742	Dilation
MB23221	MBD-FM-BB-2	MB23220	125-2744	Dilation
NT54511	NT-E-30/23-P	NT54510	133-5543	Polypectomy
PC54011	BBP-50	PC54010	129-0676	GI Accessories
PC54021	BBP-50-S	PC54020	131-2219	GI Accessories
PC54031	BBP-50-M	PC54030	131-2214	GI Accessories
PS51011	NPFS01-01023230	PS51010	129-0186	Polypectomy
PS51021	NPFS01-01523230	PS51020	129-0673	Polypectomy
PS51031	NPFS01-02023230	PS51030	129-0674	Polypectomy
PS51041	NPFS01-02423230	PS51040	129-0675	Polypectomy
PS51051	NPFS01-03023230	PS51050	129-0187	Polypectomy
PS51071	NPFS02-01023230	PS51070	129-0188	Polypectomy
PS51081	NPFS02-01523230	PS51080	129-0189	Polypectomy
PS51091	NPFS02-02023230	PS51090	129-0190	Polypectomy
PS51101	NPFS02-02423230	PS51100	129-0191	Polypectomy
PS51111	NPFS02-03023230	PS51110	129-0192	Polypectomy
PS51121	NPFS02-03623230	PS51120	129-0193	Polypectomy
PS52031	NPFS02-11523230	PS52030	129-0194	Polypectomy
PS52041	NPFS02-12523230	PS52040	129-0195	Polypectomy



Order Number	Ref #	UPN Pouch	Henry Schein Item Number	Product Family
QL22011	MDBR1A0504A	QL22010	131-2220	GI Accessories
QL22021	MDBR1B0016A	QL22020	131-2222	GI Accessories
QL22031	MDBR1B0015A	QL22030	131-2221	GI Accessories
QL22041	MDBR1B0005A	QL22040	131-2215	GI Accessories
QL22051	MDBR1B0013A	QL22050	131-2216	GI Accessories
QL22061	MDBR1B0014A	QL22060	131-2217	GI Accessories
QL6015	QL6015	QL6015	125-2750	Dilation
RC30381	ROCC-F-26-235-C-R	RC30380	128-5655	Hemostasis
RC30385	ROCC-F-26-235-C-R	RC30384	128-5657	Hemostasis
RC30411	ROCC-C-26-235-C-N	RC30410	132-5187	Hemostasis
RC30415	ROCC-C-26-235-C-N	RC30414	132-5180	Hemostasis
RC30441	ROCC-D-26-235-C-N	RC30440	132-5723	Hemostasis
RC30445	ROCC-D-26-235-C-N	RC30444	132-5724	Hemostasis
SC62241	WP-18/2200	SC62240	133-5544	Polypectomy
SC62301	WP-23/2500	SC62300	133-5545	Polypectomy

GENERAL TERMS AND CONDITIONS

Micro-Tech Endoscopy USA, Inc., dba Micro-Tech USA, Inc. ("MTUSA") sells to the Buyer subject to the following terms and conditions.

PAYMENT TERMS AND PRICE

- 1.1 All invoices submitted by MTUSA to Buyer are payable in full within 30 days after the date of such invoice. If payment is not received by MTUSA from Buyer within this 30-day period, interest shall accrue on any unpaid balance from the date of the invoice at the rate of 1.5% per month (18% per annum), but in no event at a rate greater than the maximum rate permitted by applicable law. Buyer is responsible for all collection costs, including reasonable attorney fees.
- 1.2 When pricing is offered based upon a quantity of items and Buyer does not purchase all such items as listed, MTUSA reserves the right to adjust the prices of the items actually purchased by Buyer and Buyer shall pay such adjusted prices upon demand.
- 1.3 Buyer shall pay any sales, use, excise or similar tax applicable to the purchase of the product.
- 1.4 Prices and quotations are subject to change.
- 1.5 Buyer shall be fully responsible for and will be billed by MTUSA for all freight, insurance and other transportation costs.

PURCHASE ORDERS, DELIVERY, TITLE, RISK OF LOSS, INSPECTION AND RETURNS

- 2.1 No purchase order submitted by Buyer for any product shall be binding on MTUSA unless specifically accepted by MTUSA in writing or by fulfillment thereof. Any purchase order or acknowledgment of Buyer, whether oral, printed, stamped, typed or written shall be strictly subject to and governed by these terms and conditions and none of the provisions of such purchase orders or acknowledgments shall be applicable to the extent they conflict with MTUSA's terms and conditions of sale.
- 2.2 MTUSA shall have the right, in the exercise of its sole and absolute discretion, to cancel any purchase order accepted by MTUSA or to delay shipment if Buyer fails to meet any payment schedule or supply reasonable credit or financial information requested by MTUSA.
- 2.3 MTUSA will have the right to cancel any purchase order, in whole or in part, without incurring any liability to Buyer and MTUSA shall not be liable for failure to fulfill any purchase order or to perform any contract in the event of acts beyond its control, including, without limitation, reasons of force majeure.
- 2.4 Orders are FOB shipping point. Standard shipping method will be ground delivery. Express or Overnight shipping is available at Buyer's request.
- 2.5 Delivery dates are estimates only. Buyer shall not be relieved of its obligation to purchase all of the products stated in its purchase order because of MTUSA's failure to meet any desired or requested delivery date. Title and risk of loss of any product ordered by Buyer shall pass to Buyer upon pickup and acceptance by the carrier.
- 2.6 Buyer will inspect each shipment of product and notify MTUSA within 10 days of receipt of product of any cause for rejecting the product or the product will be deemed accepted by Buyer.
- 2.7 MTUSA offers a 100% refund on all unused and unopened products. MTUSA will not accept any return of any product for credit, replacement, or refund without prior written authorization from MTUSA. Shipping and handling are not refundable.

Micro-Tech Endoscopy USA, Inc., dba Micro-Tech USA, Inc. ("MTUSA") sells to the Buyer subject to the following terms and conditions.

LIMITED WARRANTY AND DISCLAIMERS

- 3.1 MTUSA WARRANTS TO BUYER THAT, FOR THE EARLIER OF ONE YEAR FROM THE DATE OF PURCHASE OR UNTIL THE PRODUCT IS USED BY BUYER, THE PRODUCTS WILL BE FREE FROM DEFECTS IN MATERIALS AND WORKMANSHIP WHEN STORED AND USED IN ACCORDANCE WITH THE INSTRUCTIONS FOR STORAGE AND USE PROVIDED BY MTUSA AND IN ACCORDANCE WITH APPLICABLE REGULATORY REQUIREMENTS. DESCRIPTIONS OR SPECIFICATIONS APPEARING IN MTUSA'S LITERATURE ARE MEANT TO GENERALLY DESCRIBE THE PRODUCTS AND DO NOT CONSTITUTE ANY EXPRESS WARRANTIES. IN THE EVENT THAT MTUSA GIVES TECHNICAL ADVICE WITH RESPECT TO THE PRODUCT, IT IS AGREED THAT SUCH ADVICE IS GIVEN WITHOUT ANY LIABILITY ON MTUSA'S PART. ANY GUARANTEE OF SPECIFIC PROPERTIES OF THE PRODUCTS SHALL ONLY BE EFFECTIVE IF AND TO THE EXTENT SPECIFICALLY CONFIRMED BY MTUSA IN WRITING. THIS WARRANTY SHALL NOT APPLY FOR PRODUCT FAILURE OR DEFICIENCY DUE TO IMPROPER STORAGE, ALTERATION OR THE CONSEQUENCES OF USES FOR WHICH THE PRODUCTS WERE NOT DESIGNED.
- 3.2 THE WARRANTIES, OBLIGATIONS AND LIABILITIES OF MTUSA AS SET FORTH HEREIN ARE EXCLUSIVE. MTUSA DISCLAIMS AND BUYER WAIVES ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, ARISING BY LAW OR OTHERWISE, WITH RESPECT TO THE PRODUCTS AND ANY OTHER GOODS OR SERVICES DELIVERED BY MTUSA INCLUDING BUT NOT LIMITED TO ALL OTHER EXPRESS AND IMPLIED WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE; AND ANY IMPLIED WARRANTY ARISING FROM COURSE OF PERFORMANCE, COURSE OF DEALING OR USAGE OF TRADE.
- 3.3 MTUSA'S LIABILITY UNDER THIS WARRANTY IS LIMITED TO: (1) THE REPLACEMENT OF THE PRODUCTS OR THE RE-SUPPLY OF EQUIVALENT PRODUCTS; OR (2) THE REPAIR OF THE PRODUCTS OR PAYMENT OF THE COST OF REPAIRING THE PRODUCTS. MTUSA SHALL HAVE NO OBLIGATION OR LIABILITY, WHETHER ARISING IN CONTRACT (INCLUDING WARRANTY), TORT (INCLUDING ACTIVE, PASSIVE, OR IMPUTED NEGLIGENCE, STRICT LIABILITY, OR PRODUCT LIABILITY) OR OTHERWISE, FOR ANY SPECIAL, CONSEQUENTIAL, PUNITIVE, EXEMPLARY, INCIDENTAL OR INDIRECT DAMAGES.

CHOICE OF LAW AND JURISDICTION

- 4.1 All claims and disputes relating to or arising from these Terms and Conditions or the sale of MTUSA products (whether in contract, tort or otherwise and including statutory and consumer protection claims) will be governed by Michigan law without regard to conflicts of law. Buyer consents to the jurisdiction of the courts of Washtenaw County, Michigan.



MTENDOSCOPY.COM

©2018 Micro-Tech Endoscopy USA Inc. All rights reserved. All other trademarks are the property of their respective owners. MPW001 Rev 5